



**15854 LED RING LIGHT**

- 72 LED bulbs with 10 step variable intensity light control
- Sectional lighting control; allows users to direct the light from 4 separate quadrants
- 60.5mm inner diameter
- Color Temperature: 6400 Kelvin
- Lamp life: 10,000 hours
- Power supply: AC 90—265v, 50/60Hz
- Output: 12v 400mA
- RoHS compliant



**15852 LED RING LIGHT**

- 60 LED bulbs with variable light control and external power supply.
- 62mm inner diameter
- Lamp life: 10,000 hours
- Power supply: AC 90—265v, 50/60Hz
- Output: 24v 6w
- External power supply dimensions: 122 x 94 x 44 mm




## 15853 POLARIZING LED RING LIGHT

- Variable LED ring light with removable polarizer attachment.
- The polarizer reduces glare from highly reflective surfaces and is easily rotatable.
- 70 mm inner diameter
- Lamp life: 10,000 hours
- Power supply: AC 90—265v, 50/60Hz
- Output: 24v 6w
- External power supply dimensions: 122 x 94 x 44 mm



## 15856 LED RING LIGHT

- 60 LED bulbs with variable intensity light control
- ESD safe ABS plastic housing
- 60mm inner diameter
- Lamp life: 10,000 hours
- Power supply: AC 90 - 265v, 50/60Hz
- Output: 6 watt



## 15857 LED DUAL GOOSENECK ILLUMINATOR

- Dual gooseneck variable intensity LED illuminator
- Each gooseneck has a 3 watt LED bulb allowing for pinpoint light control
- Lamp life: 10,000 hours
- Power supply: AC 100-240v, 50/60Hz

# UNITRON®

73 Mall Drive, Commack, New York 11725

Tel: (631) 543-2000 Fax: (631) 589-6975

Email: [info@unitronusa.com](mailto:info@unitronusa.com) • Website: [www.unitronusa.com](http://www.unitronusa.com)

**ISO 9001  
Certification**

Design and production  
adheres to ISO9001  
international quality standard.

**ISO 14001  
Certification**

Design and production meets the  
requirements of international standard  
ISO 14001 for environmental management.