MICROSCOPES by UNITRON



LFOLED40
LED ILLUMINATOR



SPECIFICATIONS

- Engineered for superior cooling, with quiet operation
- Eliminates the continuous purchase and changing of EKE bulbs
- Half the energy (80 watts) versus the use of standard 150 watt EKE halogen F/O illuminators
- 6000K color temperature
- Dimmable brightness control
- Minimum LED life of 25000 hours
- Fits all LX accessories current or legacy
- Universal power supply for worldwide usage (96-264 VAC, 50-60Hz)
- Weight 7 lbs.
- Size: 6″ H X 4.75″ W X 9″ L
- 1 year warranty

72 Hour | QUICK SHIP

UNITRON®
73 Mall Drive, Commack, NY 11725 • 631-543-2000 (P) • 631-589-6975(F)
orders@unitronusa.com • https://microscopes.unitronusa.com